Introduction

Name: ZL-3022 Lab Hot press machine





 Tel: <u>+86-769-3365 7986</u>
 Fax: <u>+86-769-3365 7987</u>

 Website: <u>www.zltester.com</u>
 E-mail: <u>Zhongli@zltester.com</u>

Address: No.62, Xinhe Management of Area, Wanjiang District, Dongguan city, Guangdong province, China

ZhongLi Instrument Technology Co., Ltd.

Instruction:

The material is placed in the casting mold and clamped between the electric-heating flat plates to apply pressure and temperature, make the material shaping, testing and serving as the basis for the ingredient of mass production.

10~50ton select-able
Room temp. \sim 300 $^{\circ}$ C
LED Indicator,press key set input, the value indicate output
LED display,time range 0.1S \sim 99min \sim 9.9hr
Electrical Heating module
300×300mm
300×300mm
100mm
SKD chromium-molybdenum alloy
Mirror chrome plated
Room temp. $\sim +5^{\circ}\mathbb{C}$
Tap water cooling/compressor cooling water(optional)
The proportion hydraulic pressure current capacity loop circuit control
system,
the GPY oil pump driver high performance electroencephalographic
matching high accuracy non resistance piston of oil cylinder carry
out constant speed rate control mode, can p reelection one, two
times pressure and time.
Max 200mm
11mm/s
$0{\sim}350$ kg/sq.cm
10.6HP
1000×500×1500 mm
3 ∮ , AC380V, 20A

 Tel: <u>+86-769-3365 7986</u>
 Fax: <u>+86-769-3365 7987</u>

 Website: <u>www.zltester.com</u>
 E-mail: <u>Zhongli@zltester.com</u>

Address: No.62, Xinhe Management of Area, Wanjiang District, Dongguan city, Guangdong province, China



东莞中力仪器科技有限公司

ZhongLi Instrument Technology Co., Ltd.

Weight	About 658kg
--------	-------------

 Tel: <u>+86-769-3365 7986</u>
 Fax: <u>+86-769-3365 7987</u>

 Website: <u>www.zltester.com</u>
 E-mail: <u>Zhongli@zltester.com</u>

Address: No.62, Xinhe Management of Area, Wanjiang District, Dongguan city, Guangdong province, China